

Contents in Box:

1 Medical Alert System (with AC Adapter attached) - 1 Transmitter (necklace or wristband style) - 1 Telephone Cord

SETUP INSTRUCTIONS

- 1) Remove your Medical Alert System and other contents from its packaging.
- 2) Turn the Medical Alert System around so you can view the rear side of the unit. You should see two (2) ports labeled "LINE" and "PHONE".



- 3) Plug one end of the enclosed telephone cord into the port labeled "LINE".



- 4) Plug the other end of the SAME telephone cord directly into your telephone jack (usually attached to your wall). If you already have a telephone or other device plugged into the phone jack, don't worry. Simply unhook it and these instructions will allow you to reconnect it through the Medical Alert System.



NOTE: IT IS EXTREMELY IMPORTANT THAT ONE (1) TELEPHONE CORD BE CONNECTED FROM THE PORT LABELED "LINE" DIRECTLY TO YOUR TELEPHONE JACK.

- 5) Plug the AC Adapter (power cord) into any standard electrical outlet. You should hear the unit emit a "beep" when the plug is inserted into the outlet.



- 6) If you had a telephone or other device originally hooked into the phone jack, you can now plug it into the back of the Medical Alert System where the port is labeled "PHONE". Please note that any device (such as a telephone, computer, answering machine, caller ID box, etc.) can be plugged into this port.

ONCE THE SYSTEM IS CONNECTED, PLEASE TEST IT BY PRESSING THE TRANSMITTER BUTTON. IN ABOUT 40 SECONDS THE EMERGENCY OPERATOR SHOULD SPEAK TO YOU THROUGH THE SYSTEM. PLEASE TELL HER YOU ARE TESTING THE SYSTEM AND ASK HER TO READ YOUR INFORMATION (NAME, ADDRESS, PHONE NUMBER) BACK TO YOU. IF THE INFORMATION IS INCORRECT, PLEASE NOTIFY MEDSCOPE IMMEDIATELY. IF SHE DOES NOT COME ON, PLEASE CHECK YOUR CONNECTIONS OR CALL FOR ASSISTANCE.

FOR TECHNICAL ASSISTANCE OR CUSTOMER SUPPORT, CALL 1-800-645-2060